




Easy




Features & Benefits


- **Fully Rotational** - No top, bottom, front or back
- Ideal pressure relieving cushion for use in the nursing home environment.
- The Easy cushion has a layer of FloGel sandwiched between an upper and lower layer of **viscoelastic foam**
- The combination of pressure relieving **FloGel** and **viscoelastic foam** facilitates effective pressure re-distribution while allowing the patient to immerse themselves comfortably into the cushion
- Conforms to the shape of the user, increasing surface area and relieving pressure
- Clear protective polyurethane film coating is easy to wipe clean and is suitable for incontinent users
- Lightweight and easy to transport
- Includes outer cover (see page 12)

"The gel in the Easy cushion has the effect of reducing the risk of the cushion bottoming out whilst maintaining effective pressure re-distribution and facilitating the user becoming 'immersed' into the slow memory foam. Because of this design it allows the cushion to be fully rotated in use and is excellent for users and/or carers who have difficulty in applying more complex products"

International Journal of Rehabilitation, Volume 13, Issue 1, 2005



Pressure diagram with Easy cushion



Pressure diagram without wheelchair cushion

Product Code	Dimensions* (Width x Length x Height)	Cushion Weight
RE005-1416	36cm x 40cm x 5cm / 14" x 16" x 2"	1.7kg/3.8lb
RE005-1418	36cm x 45cm x 5cm / 14" x 18" x 2"	1.8kg/4.0lb
RE005-1616	40cm x 40cm x 5cm / 16" x 16" x 2"	1.8kg/4.0lb
RE005-1717	43cm x 43cm x 5cm / 17" x 17" x 2"	1.9kg/4.2lb
RE005-1816	45cm x 40cm x 5cm / 18" x 16" x 2"	1.9kg/4.2lb
RE005-1818	45cm x 45cm x 5cm / 18" x 18" x 2"	2.0kg/4.4lb
RE005-1820	45cm x 50cm x 5cm / 18" x 20" x 2"	2.1kg/4.6lb
RE005-2020	50cm x 50cm x 5cm / 20" x 20" x 2"	2.2kg/4.8lb
RE005-2218	55cm x 45cm x 5cm / 22" x 18" x 2"	2.3kg/5.0lb

The weight limit for the Easy cushion is 130kg/287lbs

*A tolerance of +/- 5mm is applied to dimensions